

# MegaCore/MC-A

## Re-usable Biopsy Gun & Needles suitable for most Biopsy Devices



### **Instructions:**

**For cocking, pull until the trigger is loaded. Pull back and rotate the actuator to retrieve the sample from the stylet.**

- High velocity of the needle for excellent histological cores
- Minimised trauma to organ tissue
- No deflection of needle and precise localisation
- Easy specimen retrieval without needle removal
- The gun is reusable, the needle is disposable

### **Reusable Gun**

Order Code	Description
B-T072222	Penetration depth 22 mm